



PCR Amplification efficiency, relative quantification, and the analysis of alterations in cellular gene expression patterns.

Olga V Petrauskene, Jennifer Maclean, Ada Wong & Manohar R Furtado, Applied Biosystems, R & D



Regulation of Gene Expression


Occurs at multiple levels. Changes in transcript levels is a key mechanism for regulation.

- Study of changes in transcript levels by relative quantification is useful if it can be done in HT mode without the need for internal controls and / or standard curves.
- Assays with near 100% amplification efficiency.
- Endogenous controls that do not change during the experiment.
- Large number of genes involved in cellular regulatory networks / pathways.
- Specificity: distinguish between transcripts from paralogous genes.




Gene Expression: 5' Nuclease Chemistry

- High sensitivity
- Large dynamic range: 10^6
- Closed tube chemistry - no post-PCR processing
- High efficiency (near 100%) associated with good designs. Useful for HT relative quantification.
- High degree of specificity from two primers and a probe.
- Contamination control is possible
- MGB enables shorter probes for greater discrimination:
 - increased fidelity to the intended target
 - increased flexibility of design




5' Nuclease Assay



1. Design Strategies in session 8.
2. Expression Testing.
3. Standards for Quantification
4. Amplification efficiency.
5. Relative Quantification (RQ).
6. Requirements for assay specificity / selectivity.
7. Application to IFN regulation.
8. Gene Silencing:RNAi & DNA methylation



Dynamic Range and Sensitivity



C_t


Sensitivity 50 –100 copies

Dynamic Range: 6-logs

confidential

Furtado Freising Talk March, 2004

Absolute Quantification: Use of Standard Curves




Wang and Morris Anal. Biochem. 269:198-201 (1999)
 Difficult to generate standard curves for large gene sets.



Co-amplification of a Competing IQC Target Identical in Size to the Target

Efficiency Tar Amp ~ 100 %

Efficiency IQC Amp ~ 100%



Introduction of 200- 500 copies/ml of IQC does not inhibit target amplification in the 50 to 5 x 10⁶ copies/mL range. Needs extensive optimization

Co amplifying targets the same size as the Target amplicon.

1. Furtado M R; Murphy R; Wolinsky S M. Quantification HIV-1 tat mRNA as a marker for assessing the efficacy of antiretroviral therapy. *Journal of Infectious Diseases*. 167(1): 213-6, 1992.
2. Mulder J, McKinney N, Christopherson C, Sninsky J, Greenfield L, Kwok S. Rapid and simple PCR assay for quantitation of human immunodeficiency virus type 1 RNA in plasma: application to acute retroviral infection. *J Clin Microbiol*. 1994;32:292-300. BASIS FOR THE ROCHE AMPLICOR ASSAY

Co amplifying targets smaller in size compared to the Target amplicon as IQCs

3. Furtado, MR, DS Callaway, JP Phair, KJ Kunstman, JL Stanton, CA Macken, AS Perelson and SM Wolinsky. Persistence of HIV-1 Transcription in Peripheral Blood Mononuclear Cells of Patients Receiving Potent Antiretroviral Therapy [see editorial review p1672]. *New England Journal of Medicine*. 340(21): 1614-1622, 1999.
4. Wang AM, Doyle MV, Mark DF. Quantitation of mRNA by the polymerase chain reaction. *Proc Natl Acad Sci U S A*. 1989;86:9717-21 . CYTOKINES

confidential

Furtado Freising Talk March, 2004

Need for Quantification Standards

Non synthetic standards

Bacteria or virus grown in culture and titered

Cellular RNA (reference RNA pools)

DNA from well characterized cell lines.

(Quantified by culture, EM , biochemical methods).

Synthetic standards

Synthetic RNA and DNA templates.

DNA cloned in plasmids.

RNA generated by *in vitro* transcription.

Recombinant retrovirus with target sequence.

Recombinant phage with target sequence.

The phage and viral based standards can be spiked into the sample at the sample preparation stage and carried through the entire process.

[NIST web site](#)

confidential

Furtado Freising Talk March, 2004

Quantification Relative to an Endogenous Control

Cellular transcripts whose levels do not change during the course of the experiment can serve as endogenous controls / reference.

18S RNA

ADA

β -actin

GAPDH

Rb

EF-1 α

UbcH5B

MLN51

Warrington et al., *Physiol. Genomics*, **2**:143-147 (2000)

Hamalainen et al., *Anal Biochem.* **299**:63-70 (2001)

Is there an endogenous control ?

User defined from experience with the system.

Needs validation.

Multiple controls.

Spike in a constant amount of an **exogenous standard RNA** to serve as an QS into all samples.

Normalize data across all samples to identify an endogenous control that changes minimally across the study.

confidential

Furtado Freising Talk March, 2004



Relative Quantification

1. Using efficiency values in the PCR equation

$$\frac{(1+E_{tar}) \{Ct_{tarcal} - Ct_{tarsam}\}}{(1+E_{ref}) \{Ct_{refcal} - Ct_{refsam}\}} = \text{Fold Change Ratio}$$

Liu and Saint, *Anal. Biochem.* **302**:52-59; 2002; Pfaffl et al., *NAR* 29:2002-2007

Efficiency values derived for target gene and reference can be used with the Ct values obtained to calculate fold change.

Accurate determination of assay efficiency is difficult. Requires the use of a broad dynamic range (5 to 6-logs) of target dilution and multiple determinations.

2. Using the $\Delta\Delta Ct$ method Livak & Schmittgen; 2001 *Methods* 25:402-8


If $E_{ref} = E_{tar} = 1$ the above equation will change to $2^{-\Delta\Delta Ct}$

(Livak 1997; ABI Prism 7700 SDS System User Bulletin # 2)

Most useful for higher throughput. No need for efficiency values.



Variability in Measuring Efficiency Values: Dependence on dilution range.



Range 2- log

82-112

97.6 Avg

N= 76

Range 5- log


92-105

98.4 Avg

N= 22



Efficiency Determination > 1000 assays




Amp Size 60-200 bp
GC content 25% to 79%
Free energy zero to -17 kcal/mole
Designed by TaqPipe; run with UMM / universal conditions

Process validation to ensure that 99.75% of the assays would have near 100% efficiency




EFF. = 94.7%

| GAPDH | | | | | | |
|-------------|-----------|----------|-------|--|---------|----------|
| DILUTION | GENE NAME | Ct | Ctdev | | | EFF % |
| 0.0001 | GAPDH | 11.235 | 8 | | -3.4568 | 94.66366 |
| 0.00001 | GAPDH | 14.22611 | 7 | | | |
| 0.000001 | GAPDH | 17.51864 | 6 | | | |
| 0.0000001 | GAPDH | 20.98957 | 5 | | | |
| 0.00000001 | GAPDH | 24.4975 | 4 | | | |
| 0.000000001 | GAPDH | 28.09658 | 3 | | | |
| 1E-10 | GAPDH | 31.57578 | 2 | | | |
| 1E-11 | GAPDH | 35.28958 | 1 | | | |




EFF. = 99.7%

| GAPDH | | | | | | |
|-------------|----------|-----------|-------|--|---------|-----------|
| DILUTION | GENE NAM | Ct | Ctdev | | | EFF % |
| 0.0001 | GAPDH | 11.234996 | 8 | | -3.3288 | 99.714021 |
| 0.00001 | GAPDH | 14.226113 | 7 | | | |
| 0.000001 | GAPDH | 17.518643 | 6 | | | |
| 0.0000001 | GAPDH | 20.989574 | 5 | | | |
| 0.00000001 | GAPDH | 24.497505 | 4 | | | |
| 0.000000001 | GAPDH | 28.096584 | 3 | | | |
| 1E-10 | GAPDH | 31.575781 | 2 | | | |
| 1E-11 | GAPDH | 35.28958 | 1 | | | |




EFF. = 97.3%

| GAPDH | | | | | | |
|-------------|----------|----------|-------|--|---------|----------|
| DILUTION | GENE NAM | Ct | Ctdev | | | EFF % |
| 0.0001 | GAPDH | 11.235 | 8 | | -3.3884 | 97.29884 |
| 0.00001 | GAPDH | 14.22611 | 7 | | | |
| 0.000001 | GAPDH | 17.51864 | 6 | | | |
| 0.0000001 | GAPDH | 20.98957 | 5 | | | |
| 0.00000001 | GAPDH | 24.4975 | 4 | | | |
| 0.000000001 | GAPDH | 28.09658 | 3 | | | |
| 1E-10 | GAPDH | 31.57578 | 2 | | | |
| 1E-11 | GAPDH | 35.28958 | 1 | | | |



EFF. = 90.6%

| GAPDH | | | | | | |
|-------------|----------|----------|-------|--|---------|----------|
| DILUTION | GENE NAM | Ct | Ctdev | | | EFF % |
| 0.0001 | GAPDH | 11.235 | 8 | | -3.5678 | 90.67105 |
| 0.00001 | GAPDH | 14.22611 | 7 | | | |
| 0.000001 | GAPDH | 17.51864 | 6 | | | |
| 0.0000001 | GAPDH | 20.98957 | 5 | | | |
| 0.00000001 | GAPDH | 24.4975 | 4 | | | |
| 0.000000001 | GAPDH | 28.09658 | 3 | | | |
| 1E-10 | GAPDH | 31.57578 | 2 | | | |
| 1E-11 | GAPDH | 35.28958 | 1 | | | |




confidential

Furtado Freising Talk March, 2004



Relative Quantification



Dilution range of human cellular RNA 100 /50/30/10/1ng Spike in a constant amount of an external transcript into all tubes in pool dilution series. Non-human transcript. Determine FC & compare with expected FC.



RQ Experimental Design; Reproducibility

1. Number of Assays Selected 15
 - Low Expressors 26-30 Ct 10,000 to 2000 5
 - Moderate Expressors 22-26 Ct 100,000 to 10,000 5
 - High Expressors 18-22 Ct 1,000,000 to 100,000 5
2. Dilutions/ Amounts: 100; 50; 30; 10 and 1 ng /well
Spike in a constant amount of external transcript
3. Number of Operators 3
4. Number of Runs 3
5. Number of sites / instrument 1
6. Number of Replicates 3



RQ: Results Summary

Operator to Operator Variation: 2-Fold Change

| | N | Range | AVG | SD | Median | 95% CI | < 1.5 | Accuracy (1.33) | Accuracy (1.50) |
|--------------------------|-----|-----------------|------|------|--------|-----------------|-----------|--------------------|--------------------|
| Exp_FC 2 and All OPER | 387 | 0.97 to 6.03 | 2.23 | 0.5 | 2.18 | 1.23 to 3.23 | 6 or 1.5% | 97.40% | 96.80% |
| Exp_FC 2 and OPER 1 | 120 | 1.59 to 4.71 | 2.29 | 0.53 | 2.15 | 1.23 to 3.35 | 0 | 100% | 100% |
| Exp_FC 2 and OPER 2 | 133 | 0.97 to 6.03 | 2.14 | 0.62 | 2.03 | 0.90 to 3.38 | 6 or 4.5% | 92.50% | 90.90% |
| Exp_FC 2 and OPER3 | 134 | 1.76 to 3.14 | 2.29 | 0.29 | 2.23 | 1.71 to 2.87 | 0 | 100% | 100% |

Robust RQ results for FC in the 1.66 to 100 range.

RQ with the $\Delta\Delta C_t$ method is more accurate than methods using efficiency values in the PCR equation.

confidential

Furtado Freising Talk March, 2004



Summary : Quantification


- Internal calibrators or standard curves are necessary for accurate absolute quantification.
- For relative quantification it is useful to have an endogenous control that is minimally altered during the study / experiment.
- Quantification standards that can be added into the sample (or during NA extraction) and carried through the entire process are useful for normalization across large studies and identification of appropriate endogenous controls.



Conclusions .. continued

- Efficiency measurements have a high degree of variability. Using these values in RQ measurements gives less robust results.
- The $\Delta\Delta C_t$ method produced a better correlation to expected results as compared to measurements using efficiency values in PCR equations.
- Good assay design strategies can ensure near 100% efficiency for real time PCR assays.
- Relative quantification using the $\Delta\Delta C_t$ method generated reliable results and is amenable to HT processing.

Assay Specificity CYP2A13 vs CYP 2A6



Selectivity by Primer Mismatches


B

| | Exon 8 | | | Exon 9 | |
|----------------|------------|------------|------------|------------|----------------|
| <i>CYP2D6</i> | GCCTTCCTGC | CTTTCTCAGC | AGGCCGCCGT | GCATGCCTCG | GGG.AGCCC |
| <i>CYP2D7P</i> | GCCTTCCTGC | CTTTCTCAGC | AGGCCGCCGT | GCATGCCTCG | GGG.AGCCC |
| <i>CYP2D8P</i> | GCCTTCCTGC | CTTTCTCAGC | AGGCCGCCGT | GCATGCCTCG | GGCCAGCCC * |
| <hr/> | | | | | |
| <i>CYP2D6</i> | CTGGCCCGCA | TGGAGCTCTT | CCTCTTCTTC | ACCTCCCTGC | TGCAGCACI |
| <i>CYP2D7P</i> | CTGGCCCGCA | TGGAGCTCTT | CCTCTTCTTC | ACCTCCCTGC | TGCAGCACI |
| <i>CYP2D8P</i> | CTGGCCCGCA | TAGAGCTCTT | CCTCTTCTTC | ACCTCCCTGC | TGCAGCACI * |
| <hr/> | | | | | |
| <i>CYP2D6</i> | TCAGCTTCTC | GGTGCCCACT | GGACAGCCC | GGCCCAGCCA | CCATGGTGT |
| <i>CYP2D7P</i> | TCAGCTTCTC | CGTGGCCGCC | GGACAGCCC | GGCCCAGCCA | CTCTCGTGT |
| <i>CYP2D8P</i> | TCAGCTTCTC | GGTGCCCACT | GGACAGCCC | GGCCCAGCCA | CTCTCGTGT |
| | | * * * * | | | * * * |
| <hr/> | | | | | |
| <i>CYP2D6</i> | CTTTGCTTTC | CTGGTGAGCC | CATCCCCCTA | TGAGCTTTGT | GCTGTGCC |
| <i>CYP2D7P</i> | CGTCAGCTTT | CTGGTGAGCC | CATCCCCCTA | CGAGCTTTGT | GCTGTGCC |
| <i>CYP2D8P</i> | CGTCGGCTTT | CTGGTGAGCC | CATCCCCCTA | TGAGCTTTGT | GCTGTGCC |
| | * * * * * | ** | | * | |

Anal Biochem. 300:
121-131. 2002
Cytochrome P450
CYP2D6



Effect of a Mismatch on Ct



The Ct value is dependent on both PCR amplification & probe cleavage.

Effect of mismatches in primers & probe on Ct.

Position specific effects.

Mismatch specific effects.


Combinatorial effects of multiple mismatches

Relative contributions from primers & probes.

Over 3,000 designs tested to generate data.



In silico QC against the known transcriptome.

Selective TaqPipe: Bioinformatics Tool



Testing a large number of designs to define contributions to selectivity by **mismatch on a primer and a probe: type and location**

Determining Level of Assay Specificity





Conclusions

➤ Real time PCR assays that target transcripts emanating from closely related genes should have high selectivity.


➤ Suppression of co-amplification can be achieved by placing primers at positions where there are mismatches.

➤ Mismatches within probe binding regions will also enhance selectivity.


➤ The level of selectivity should be defined by the user based on the study / experiment and is related to transcript abundance and levels of induction / suppression of the related genes.



IFN Pathway



Comparison of Interferon-Induced Fold Change in Five Genes





| Gene Name | Fold Change at 24 hours (Literature) | Fold Change at 24 hours in study |
|------------------|---|---|
| PRKR | 7 | 7.03 |
| L22342 | 10 | 12.87 |
| IFI16 | 12.7 | 15.1 |
| SSA1 | 6.6 | 18.8 |
| ISGF3/IRF9 | 9.6 | 12.7 |
| STAT1 | 11.2 | 22.3 |
| OAS1 | 25.2 | 29.6 |
| IFIT1 | 124.6 | 490.12 |

Yuen et al., NAR 30 No 10 e48 ; 2002

Arrays

PCR


Tested over a 100 genes with good correlation to literature.

Microarray analysis, mass spec analysis for protein profile and PTMs in progress



BRCA1 (3 siRNA vendors):

AoD (3') and AbD (5') assays detect 30-85% RNA knock-down by different siRNAs




HeLa cells
Control


AbD (exon 2-3)

AoD (exon 15-16)


Exon:
NM_007294
(7100 bps)



HeLa cells
siRNA treated
24 Hours




BRCA1 gene: 23 splice isoforms




26 exons and 2 genes found in this genomic region spanning 81891 bp.

[View graphic only](#)



CyclinD1: siRNA



More than 90% knock-down detected



Acknowledgements

Karl Guegler

Ken Livak

Dennis Gilbert

Gene Spier

Olga Petrauskene

Kathleen Shelton

Karen Poulter

Jennifer MacLean

Rixun Fang

Jeff Peralta

Michael Malicdem

Ben Curzon

Ada Wong

Sonny Young

Eileen Westphale



Labeling

For Research Use Only. Not for use in diagnostic procedures.

The PCR process and 5' nuclease process are covered by patents owned by Roche Molecular Systems, Inc. and F. Hoffmann-La Roche Ltd.

Applied Biosystems is a registered trademark and AB (Design), Applera, Assays-by-Design, Assays-on-Demand, and ViroSeq are trademarks of Applera Corporation or its subsidiaries in the US and/or certain other countries.

© 2003 Applied Biosystems. All rights reserved.